






<b>G-Ceram ZF Base - Dentine</b>										<b>A0 - D4</b>		<ul style="list-style-type: none"> <li>▶ Increase the opacity effect of dentine.</li> <li>▶ Mainshade material to build up the ceramic restoration</li> </ul>						
A0	A1	A2	A3	A3,5	A4	B0	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
<b>G-Ceram ZF Dentine</b>										<b>A0 - D4</b>		<ul style="list-style-type: none"> <li>▶ Used to increase the chroma effect of dentine layer</li> <li>▶ Mainshade material to build up the ceramic restoration according to natural teeth</li> </ul>						
A0	A1	A2	A3	A3,5	A4	B0	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
<b>G-Ceram ZF Incisal</b>										<ul style="list-style-type: none"> <li>▶ Last layer of build up the ceramic restoration according to natural teeth</li> </ul>								
I1	I2	I3	I4	I5	I6	Opal incisal clear		Opal incisal light		Opal incisal medium		Transparent						

<b>G-Ceram ZF Shoulder</b>		Bleach	<ul style="list-style-type: none"> <li>▶ Used for optimizing the marginal areas of zirconia copings.</li> </ul>
		Light	
		Medium-Light	
		Medium	
		Medium-Dark	
		Dark	












<p><b>G-Ceram ZF Dentine Modifier</b></p>		White	<p>▶ Used for preparing mamelons and creating special characterization</p> <p>▶ Used for coloring the inside to obtain effects of depth</p> <p>▶ Used for optimizing intensive dentine shade</p>
		Yellow	
		Ocher	
		Brown	
		Gray	
		Pink	
		Orange	
		Blue	
		Violet	




<p><b>G-Ceram ZF Gum Shade</b></p>		Light	<p>▶ Gum materials for esthetic appearance</p>
		Medium	
		Dark	

<b>G-Ceram ZF Incisal Modifier</b>		White	<ul style="list-style-type: none"> <li>▶ For different characterization and individualization</li> <li>▶ Provide different colors of translucency effects</li> </ul>
		Yellow	
		Ocher	
		Brown	
		Gray	
		Pink	
		Orange	
		Blue	
		Violet	

<b>G-Ceram Liquids</b>		Modelling	▶ Used to mix all dentine, incisal and additional materials
		Glaze & Stain	▶ Used to mix all glaze & stain powders

<b>G-Ceram ZF Low Glaze</b>		▶ For polish surface of restoration
-----------------------------	--	-------------------------------------

<p>G-Ceram ZF Stain</p>		White	<p>▶ Used to provide the same characteristic image as natural teeth</p> <p>▶ Used for restoration make up alone or mixed with one another</p>
		Yellow	
		Ocher	
		Brown	
		Gray	
		Pink	
		Orange	
		Blue	
		Violet	
		Green	
		Black	

<p>G-Ceram ZF Add On</p>		Light	<p>▶ Used for corrections and adjustment of shape after glaze firing</p>
		Medium	
		Dark	